

DISCOVER NON-AMBIENT X-RAY DIFFRACTION



SPEAKER: ELAA BEN FREDJ – PRODUCT SPECIALIST

JUNE 14, 2026

CANADIAN POWDER DIFFRACTION WORKSHOP 19 (CPDW 19)


AGENDA

- › Basic principles of Non-Ambient X-ray diffraction
- › The typical set-up
- › Practical considerations in NA-XRD
 - › Temperature accuracy and validation
 - › Sample height variations
- › Application overview



BASIC PRINCIPLES OF NA-XRD

- Non-ambient X-ray diffraction (NA-XRD) is an indispensable technique to understand the behavior of materials under non-ambient conditions
- Material properties may change significantly when exposed to non-ambient conditions.




Influence on:

- Crystalline Structure
- Composition
- Crystallite Size
- Internal Stresses
- Surface Layer properties


BASIC PRINCIPLES OF NA-XRD

NA-XRD allows to study a variety of material changes in-situ.

Processes	Application fields
Material formation and structures	Alloys, drug API's, catalysts, minerals..
Annealing due to heat treatment	Alloys, ceramics, polymers..
Calcination and sintering	Catalysts, MOF's, zeolites..
(De-) hydration processes	Building materials, pharmaceuticals, food industry..
Material changes during operation	Catalysts, refractory materials, batteries..



TYPICAL SET-UP FOR NA-XRD




NA-XRD – HEATER DESIGNS

Environmental Heaters

- Heat Transfer from all sides, via surrounding heater
- Capillary, spinner
- RT up to 1500°C

SAMPLE placement



Direct Heaters

- Heat Transfer from One side, with heater below the sample

Strip Heaters:

- RT up to 2300°C

Heating Plates / Heater-Cooler Combi

- - 190°C up to 1100°C

SAMPLE placement




NA XRD – HEATER DESIGNS

Environmental Heaters	Direct Heaters
<u>Advantages</u>	
✓ Good temperature homogeneity	✓ Very high temperatures possible
✓ More accurate sample temperature measurement	✓ Very fast heating / cooling rates
✓ Easy sample exchange	✓ Simpler design
✓ Inert sample holders, Sample spinner available	✓ Temperature-control less dependent on the atmosphere in chamber
<u>Challenges</u>	
✗ Thermal expansion	✗ Sample preparation difficult, reactions between sample and sample holder possible (especially for strip heater)
✗ Slower heating / cooling rates	✗ Less uniform sample temperatures
✗ Lower upper T-limit	✗ Larger sample-temperature deviations


NA XRD – TEMPERATURE MEASUREMENT

TEMPERATURE SENSOR (RTDS OR THERMOCOUPLES) POSITION



Close to the sample but without contact:

- Sample holder rotation possible




Inside / On the sample holder:

- Thermal contact between sensor and sample holder
- Heat transfer from sample holder to sample

NA-XRD – Measurement Strategies


- Temperature-dependent measurements

e.g. phase-transition



- Time-dependent measurements


e.g. reaction kinetics



PRACTICAL CONSIDERATIONS IN NA-XRD


SAMPLE TEMPERATURE ACCURACY

- Specified temperatures always refer to the temperature of the sample holder and **NOT** the sample
- The location of the temperature measurement and the area under investigation are different
- Deviation between the measured temperature and the sample surface temperature is expected





AREA UNDER INVESTIGATION

TEMPERATURE MEASUREMENT



SAMPLE TEMPERATURE ACCURACY INFLUENCES ON TEMPERATURE DEVIATION

- Design of non-ambient XRD attachment
 - Direct or environmental heater
- Applied atmosphere
 - Vacuum or gas atmosphere
- Thermal properties of the sample
 - Heat conductivity, sample thickness
- Thermal contact
 - Surface between heater and sample is crucial




NA XRD – TEMPERATURE VALIDATION

- Temperature validation: Correlating the sensor temperature with the true sample surface temperature

Validation methods:

1. Phase transition: Monitor a crystallographic transformation or sample melting and compare the measured transition temperature with tabulated values

2. Thermal expansion: Measure the change of the lattice parameter(s) during heating / cooling and compare it with thermal expansion curves from literature




NA XRD – TEMPERATURE VALIDATION

- Validation must be done under the same conditions as the measurements (gas type, pressure, sample thickness, ...)
- The validation/ reference sample should have similar thermal properties as the sample (color, type of material, ...)
- Validation only makes sense if the temperature deviation is reproducible.

Substance	Transition Temp [°C]:	Type ¹	Lit. ²	Substance	Transition Temp [°C]:	Type ¹	Lit. ²
NH ₄ NO ₃	54	p.t.	1	SiO ₂	573	p.t.	1
NH ₄ NO ₃	127	p.t.	1	Li ₂ SO ₄	577,85	p.t.	3
KNO ₃	128,7	p.t.	2	K ₂ SO ₄	584	p.t.	1
In	156,5985	m.p.	3	Sb	630,5	m.p.	4
RbNO ₃	166	p.t.	2	Rb ₂ SO ₄	653	p.t.	5
RbNO ₃	222,7	p.t.	2	Al	660,323	m.p.	3
Sn	231,928	m.p.	3	KCl	776	m.p.	4
Bi	271,4	m.p.	2	NaCl	804	m.p.	4
RbNO ₃	285	p.t.	1	Bi ₂ O ₃	820	m.p.	4
KClO ₄	299,4	p.t.	2	Ag	961,78	m.p.	3
Cd	321	m.p.	2	NaF	988	m.p.	4
Pb	327,46	m.p.	3	Au	1064,18	m.p.	3
KNO ₃	334	m.p.	2	K ₂ SO ₄	1069	m.p.	4
Zn	419,527	m.p.	3	Cu	1083	m.p.	4
AgSO ₄	426,4	p.t.	2	CaF ₂	1360	m.p.	4
CuCl	430	m.p.	4	Ca ₂ SiO ₄	1425	p.t.	4
CsCl	476	p.t.	1	Fe	1535	m.p.	4

Validation methods:

Phase transition – Well known references



NA XRD – TEMPERATURE VALIDATION

Phase Transition

Thermal Expansion

Advantages

✓ No special requirements with regards to data quality

✓ Straight-forward data evaluation

✓ Complete temperature range can be validated with one sample material

✓ Method can be applied up to very high temperatures

Challenges

✗ Only one (max. a few) temperature value per reference material

✗ Limited availability of feasible phase transformations (especially at low temperatures)

✗ Data evaluation can be complex

✗ Inaccurate at low temperatures (< 400°C), due to small lattice expansion

THERMAL EXPANSION OF THE SAMPLE HOLDER

Effect of sample height variation:


- Heating leads to thermal expansion of materials
- Thermal expansion of sample holder results in misalignment of the sample height → peak shift in diffraction pattern

$$\Delta 2\theta = \frac{-2s \cos \theta}{R}$$

s... magnitude of displacement

R... radius of goniometer

θ ... scattering angle in radians



SAMPLE HEIGHT VARIATIONS



Fig.: Typical sample holder for HTK 1200N

Possible approaches:

- ! Ignore the effect
- ! Correct during data analysis
- ! Use a parallel beam geometry

Recommended:

- ✓ **Correct the displacement during the measurement**
 - › Manual re-alignment
 - › Automated Z-stage compensation

SAMPLE HEIGHT VARIATIONS REALIGNMENT (Z-ALIGNMENT)

Correct height alignment is highly important for good diffraction data!

- This is done by lowering the complete NA-XRD attachment (importance of height-alignment stage) according to a predefined correlation table (instrument specific) or temperature vs. height change.




Fig.: Typical sample holder for HTK 1200N

APPLICATIONS

STRUCTURAL CHANGES IN BENTONITE DURING HEATING

APPLICATION 1


- Bentonite is a naturally occurring swelling clay that can absorb large amounts of water.
- Its application include drilling muds, casting sands, ceramics, catalysts, adsorbents and geological barrier materials.
- In-situ XRD measurements with an environmental heater monitors structural and phase changes during heating up to 1000 °C.




HTK 1200N: environmental heater with high temperature homogeneity _ Temperature range: 25 °C to 1200 °C
Atmospheres: Air, inert gas, vacuum (10^{-4} mbar), O₂

STRUCTURAL CHANGES IN BENTONITE DURING HEATING

APPLICATION 1



Diffractogram of the starting material consisting of montmorillonite, quartz and muscovite.



Selected diffractograms from the in-situ experiment.



EFFECT OF HYDROPHOBIC SURFACES ON WATER FREEZING

APPLICATION 2

- › Water covers approx. 70% of the planet and a significant amount of water exists as vapor in the atmosphere.
- › Water is the only non-metallic substance that expands as it freezes.
- › Expanding water can damage materials or surfaces and can lead to bursting of containers or pipes
- › Frozen roads also lead to many car accidents every year.




EFFECT OF HYDROPHOBIC SURFACES ON WATER FREEZING

APPLICATION 2

Experiment

- › The correlation between hydrophobicity of surfaces and the freezing point of water droplets deposited on the surface was investigated:
- › Experiments were performed with low temperature chamber with liquid nitrogen cooling from -5 °C to -25 °C.




TTK 600: Low-Temperature Chamber
Temperature range: -190°C to 600°C.
Measurements can be performed in air, inert gas
or high vacuum (10⁻⁴ mbar)

EFFECT OF HYDROPHOBIC SURFACES ON WATER FREEZING

APPLICATION 2


- › Silicon (Si) wafers were used as substrates.
 - › One plate was untreated and served as reference to determine the change of the freezing point.
 - › On silicon plate featured hydrophobic surface coating (grad100)
- › Water droplet (0.1 mL) was deposited on surface.



EFFECT OF HYDROPHOBIC SURFACES ON WATER FREEZING

APPLICATION 2


- › As expected, no crystalline reflections of water are visible at 20 °C.
- › The freezing point of water on the **untreated surface** is at -9 °C.
- › The reflection labeled * is most likely Si 200, which is often observed in (111) oriented Si wafers.



EFFECT OF HYDROPHOBIC SURFACES ON WATER FREEZING


APPLICATION 2

Untreated silicon




The freezing point of water on the **untreated surface** is at **-9 °C**.

Hydrophobic silicon



The freezing point of water on the **hydrophobic surface** is at **-22 °C**.




OPERANDO XRD UNDER NON-AMBIENT CONDITIONS

APPLICATION 3

Objective:

Monitor XRD peak shifts during electrochemical charge cycling.

Evaluate structural changes over a temperature range of $-5\text{ }^{\circ}\text{C}$ to $55\text{ }^{\circ}\text{C}$.



Battery sample holder (Reflection)


Measurement configuration:

- Cu X-ray tube
- Bragg-Brentano geometry using a monochromator
- External potentiostat for for operando charge/discharge cycling

Sample:


LiCoO₂ (LCO) half-cell in coin cell format

- Li metal anode
- Positive electrode: aluminum foil coated with LCO



TTK 600 Low Temperature Chamber for in situ and in operando XRD studies with temperature range of -180 to $+130\text{ }^{\circ}\text{C}$.

PEAK SHIFTS DURING BATTERY CYCLING




Heat maps of XRD scans for 003, 101, and 104 reflection of LCO during battery cycling (left), and charging/discharging curve (right).

In situ XRD monitor the structural changes in real time:

- Peak shifts
- Lattice parameter changes
- Phase transitions



NON-AMBIENT CYCLING OF LCO

- Capacity decreases at low temperatures
- Capacity is relatively stable between 10 °C and 45 °C
- At higher temperatures:
 - capacity measured during discharge drops
 - capacity obtained during charging seems to increase.



Capacity and Coulomb efficiency of LCO coin cell at different temperatures.

NA CYCLING OF LCO



Temperature-dependent contour plots of XRD scans for 003, 101 and 104 reflection of LCO during battery cycling.

NA XRD – summary

- Non-ambient XRD enables the investigation of structural changes under realistic operating conditions.
- Reliable measurements require careful control of temperature, sample alignment, and experimental conditions.
- From heating and cooling to operando battery cycling, non-ambient XRD provides direct insight into the relationship between structure and material performance.



NON-AMBIENT XRD SPECIFICATION RANGE

COVERING A HUGE RANGE OF NON-AMBIENT CONDITIONS

