

# Operando diffraction mapping of commercial Li-ion batteries

*Doing XRD with complex samples, difficult data, and spatial/temporal changes*

Canadian Powder Diffraction Workshop – July 31, 2025

**Toby Bond**

Senior Scientist: Energy Storage and X-ray Imaging  
Industry Services Group  
Canadian Light Source



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)

## Background

- Li-ion batteries: electrodes and cells
- Degradation and microcracking
- In-situ vs commercial cells

## Time-resolved mapping with Synchrotron XRD

- Experimental considerations
- Data analysis





Canadian  
Light  
Source    Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY

# Li-ion batteries: Overview


- Most common negative electrode is graphite
- Many different positive electrodes – here we only consider NMC (Nickel Manganese Cobalt oxide) in particular,  $\text{LiNi}_{0.6}\text{Mn}_{0.2}\text{Co}_{0.2}\text{O}_2$  (NMC622)





# Li-ion batteries: packs, cells, and electrodes


## Battery packs and modules



Cylindrical cell

Pouch cell

Prismatic cell



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron




**DALHOUSIE**  
UNIVERSITY

Image sources: <https://www.flickr.com/photos/argonne/14547144640>, <https://iopscience.iop.org/article/10.1149/2.0691916jes>


THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](http://lightsource.ca)

# Li-ion batteries: cells and electrode microstructure

## Graphite anode



## NMC cathode



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY



Image sources: <https://www.flickr.com/photos/argonne/14547144640>, Jeschull et al., J. Electrochem. Soc., 167, 100535 (2020)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](http://lightsource.ca)




# Cathode microcracking

- Anisotropic volume change in NMC and NCA causes significant stress at grain boundaries




# Cathode microcracking: electrical disconnection

Conductive-binder particles




Disconnected particles



# In-situ cells vs commercial cells

## In-situ cells



### Pros:

- Optimized for x-ray experiments
- can be made in house
- Well suited to short, simple experiments

### Cons:

- Not fully representative of real-world conditions
- Not built for long-term cycling
- Sensitive to alignment
- Can have poor reproducibility

## Commercial cells



### Pros:

- Real-world form factor
- Built for long-term cycling
- Highly reproducible
- Can test cells from the field

### Cons:

- Not optimized for x-ray experiments
- Non-active components can interfere with data collection
- Need high x-ray energy to penetrate cells (especially steel casings)



Canadian Light Source  
Centre canadien de rayonnement synchrotron




DALHOUSIE  
UNIVERSITY


Figure by: A Llewellyn et al., Condens. Matter, 5(4), 75 (2020)



# Long-Term Cycling of NMC622/NG Pouch Cells




- **NMC622** battery cells were cycled from 0% to 100% for 2.5 years
  - **Heavily cycled cell: deep cycling (0%-100%) for 2400 cycles (960k km) over 2.5 years**
  - **Lightly cycled cell: shallow cycling (75%-100%) for 3800 cycles (380k kms) over 1.5 years**
- A “fresh” cell was also used as a control






# Microcracking after heavy cycling


Control cell



Lightly cycled cell



Heavily cycled cell



Canadian  
Light  
Source    Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY

[T. Bond et al., J. Electrochem. Soc. 169, 080531 \(2022\)](#)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)



## Background

- Li-ion batteries: electrodes and cells
- Degradation and microcracking
- In-situ vs commercial cells

## Time-resolved mapping with Synchrotron XRD



- Experimental considerations
- Data analysis





# Synchrotron XRD of cycled pouch cells

- Flat pouch cells are great for transmission XRD
- This work was done on the BXDS high-energy wiggler beamline
- Beam width  $\sim 200\text{ }\mu\text{m}$



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron




**DALHOUSIE**  
UNIVERSITY

[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)


THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](http://lightsource.ca)

# NMC (113) and graphite (001) peaks

## NMC (113) peak






## Graphite (001) peak




# Thick samples require large sample-detector distance

- We want to reduce peak broadening due to sample thickness
- You can do this by increasing sample-detector distance, but you also need to increase beam energy (reduce wavelength)
- We looked at two configurations:
  - Detector distance = 570 mm, Beam Energy = 35 keV
  - Detector distance = 1400 mm, Beam Energy = 55 keV





# Near-equilibrium: XRD mapping after full charge/discharge



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron





**DALHOUSIE**  
UNIVERSITY

[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)

# Mapping the distribution of inactive cathode material



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY



[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)




# Non-equilibrium: operando XRD mapping

- If we repeat this raster scanning, we get both spatially and temporally resolved data
- Gray lines denote peaks that appear to be static






# Peak fitting of time-resolved XRD data



# Spatial distribution of model components at equilibrium



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY


[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)




# Spatial distribution of model components at equilibrium

## Peak fitting method




## Peak overlap method





# Time-resolved XRD mapping during Open-circuit relaxation



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron




**DALHOUSIE**  
UNIVERSITY

T. Bond et al., J. Electrochem. Soc. 171, 110514 (2024)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)

# Open-circuit relaxation: control cell



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron




**DALHOUSIE**  
UNIVERSITY

[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)



# Open-circuit relaxation: heavily cycled cell



Canadian  
Light  
Source

Centre canadien  
de rayonnement  
synchrotron





**DALHOUSIE**  
UNIVERSITY

[T. Bond et al., J. Electrochem. Soc. 171, 110514 \(2024\)](#)




THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)

# OCV relaxation: comparison





# OCV relaxation: 2D mapping



Pearson  
Correlation  
Coefficient:

$$\rho_{X,Y} = \frac{\text{cov}(X, Y)}{\sigma_X \sigma_Y}$$



Canadian  
Light  
Source  
Centre canadien  
de rayonnement  
synchrotron



DALHOUSIE  
UNIVERSITY

T. Bond et al., J. Electrochem. Soc. 171, 110514 (2024)

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](https://lightsource.ca)

# Summary

- Real-world materials and devices can be very complicated
- You're often working with incomplete or messy data that requires flexibility
- You need time to experiment with both data collection and analysis
- Synchrotron XRD is a great tool for this kind of work



Canadian  
Light  
Source    Centre canadien  
de rayonnement  
synchrotron



**DALHOUSIE**  
UNIVERSITY

THE BRIGHTEST LIGHT IN CANADA | [lightsource.ca](http://lightsource.ca)



# Thank you

## Contributing authors:

- Jeff Dahn, Dalhousie
- Sergei Gasilov, CLS
- Graham King, CLS
- Roby Gauthier, Dalhousie
- Ahmed Eldesoky, Dalhousie
- Reid Dressler, Dalhousie
- Jeffin James Abraham, Dalhousie
- Jessie Harlow, Dalhousie

This work is supported by NSERC and Tesla Canada under the Alliance grant program

The Canadian Light Source is supported by the Canada Foundation for Innovation (CFI), NSERC, the National Research Council (NRC), the Canadian Institutes of Health Research (CIHR), the Government of Saskatchewan, and the University of Saskatchewan

